

# Urine Pre-Clinical CartiLaps® ELISA

FOR RESEARCH USE ONLY

## Cartilage Degradation Marker (CTX-II) For Preclinical Research

### Animal models in arthritis

- Analyzing cartilage turnover
- Monitoring chondroprotective effects of drugs

### Screening of chondroprotective drugs

- Objective assessment of *in vitro* cartilage degradation

The Urine Pre-Clinical CartiLaps® ELISA is used for quantitative assessment of cartilage degradation. The assay detects C-telopeptide fragments of collagen type II (CTX-II) generated during cartilage erosion in animal models and cartilage cultures.

X-ray by The Royal Veterinary and Agricultural University, Dept. of Clinical Studies

nordicbioscience  
diagnostics

[www.nbdiagnostics.com](http://www.nbdiagnostics.com)

Urine Pre-Clinical CartiLaps®  
ELISA

N-MID® Osteocalcin  
ELISA

ALPHA CrossLaps®  
ELISA

Rat-MID™ Osteocalcin  
ELISA

Urine CartiLaps®  
ELISA

RatLaps™  
ELISA

Serum Pre-Clinical CartiLaps®  
ELISA


Serum CrossLaps®  
ELISA

Urine BETA CrossLaps®  
ELISA

CrossLaps® for Culture  
ELISA

# Urine Pre-Clinical CartiLaps® ELISA

ENZYME IMMUNOASSAY FOR QUANTITATIVE ASSESSMENT OF CARTILAGE DEGRADATION



## Performance Characteristics

- Method: • Competitive ELISA
- Format: • 96-well microplate with reagents sufficient to test 40 samples in duplicate
- Detection limit: • 0.75 µg/L
- Analyte: • A sequence (EKGPDP) specific for a part of the C-terminal telopeptide  $\alpha 1$  chain of type II collagen (CTX-II)
- Specimen: • Urine
- Specimen volume: • 10 µl
- Specimen CV intraassay: • <5%
- Specimen CV interassay: • <9%
- Species reactivity: • Non-human primates, bovine, horses, pigs, rabbits, rats and mice

## Sampling

Collect either spot urine or 24 hours urine. Freeze samples for prolonged storage (<math> <18^{\circ}\text{C}</math>). Do not add HCl to urine samples.

The Urine Pre-Clinical CartiLaps® ELISA kit is for in vitro use only.  
Product number #2CAL4000

## Assay Procedure

- Prior to use, prepare and equilibrate all solutions to room temperature.
- Pre-dilute urine samples 1+3 in Standard A.
- Pre-incubate by adding 100 µL **Biotinylated Antigen** to each well, cover with sealing tape, and incubate for  $30 \pm 5$  minutes at room temperature ( $18-22^{\circ}\text{C}$ ).
- Wash the immunostrips 5 times manually with **Washing Buffer** diluted 1+50 in distilled water.
- Pipette 10 µL of either **Standards**, **Control**, or unknown samples into appropriate wells followed by 150 µL **Primary Antibody**. Cover the immunostrips with sealing tape and incubate for  $21 \pm 3$  hours at  $2-8^{\circ}\text{C}$ .
- Wash.
- Add 100 µL **Peroxidase-Conjugated Antibody** solution to each well. Cover with immunostrips with sealing tape and incubate for  $60 \pm 5$  minutes in the dark at room temperature ( $18-22^{\circ}\text{C}$ ).
- Wash
- Pipette 100 µL of the **Substrate Solution** into each well and incubate for  $15 \pm 2$  minutes at room temperature ( $18-22^{\circ}\text{C}$ ) in the dark. Use sealing tape.
- Pipette 100 µL of the **Stopping Solution** into each well.
- Measure the absorbance at 450 nm with 650 nm as reference within two hours.

Species	Mean (95% CI)
Dog (study 1)	69.5 mg/mmolCr (4.1-260.4)
Dog (study 2)	666 mg/mmolCr (26.6-2132)
Guinea Pig (study 3)	43.9 mg/mmol Cr (1.1-231.2)
Guinea Pig (study 4)	388 mg/mmol Cr (57.4-840)
Rabbit (study 5)	41.9 mg/mmol Cr (7.6-114)
Rabbit (study 6)	625.2 mg/mmol Cr (60.8-1745)
Rat (study 7)	24.3 mg/mmol Cr (7.8-68.2)
Rat (study 8)	95.4 mg/mmol Cr (35.5-135)
Cell Culture Supernatant	69.8 mg/mmol Cr (2.0-237)

LITERATURE: 1. CHRISTGAU, S. ET AL. MENOPAUSE (2004); 11:508-518. 2. CHRISTGAU, S. ET AL. BONE (2001); 29: 209-215. 3. HOEGH-ANDERSEN, P. ET AL. ANNALS OF RHEUM DIS. (2004); 6(2): R169-80. 4. REIJMAN, M. ET AL. ARTHRITIS & RHEUM. (2004); 50:2471-2478. 5. ROY-BEAUDRY, M. ARTHRITIS & RHEUM (2003); 48:2855-2864.

all the way

FROM RESEARCH TO PATIENT MONITORING

Nordic Bioscience Diagnostics A/S • Herlev Hovedgade 207 • 2730 Herlev • Denmark • www.nbdiagnostics.com

ISO 9001 certified